



Exa® PACS/RIS

1.4.33_P4

Customer Release Notes

Table of Contents

New Features	3
Customer Requested Improvements	3
Other Improvements	3
Bug Fixes	3

New Features

Summary

- No new features in this release.

Customer Requested Improvements

Summary

- When locally caching studies, caching now occurs in order based on Stat level: higher Stat levels to lower Stat levels, followed by no Stat levels.
- OpalForwarder now includes more detailed logging to track files from receipt and processing to sending.
- Logs now include more information about handling of Opal packaged files by Opal Listener.
- In OpalForwarder, you can now specify a forwarding destination to use in cases when no routing match is found.

Other Improvements

Summary

- Various improvements to performance and the UI.

Bug Fixes

Summary

- Worklist filters now display all studies that match the filtering criteria, regardless of whether they are locked by another user, or whether they are DICOM studies.
- When viewing a schedule book filter, after leaving the schedule book and then returning, the default view of the schedule book now appears rather than the filter.
- When uninstalling and then reinstalling local cache, the prefetcher no longer runs before local cache is activated.
- When two Exa platform servers contain different studies with the same DICOM study ID, and a radiologist uploads a dictation for one of those studies, the dictation now uploads to the study on the correct Exa platform server.
- When the default insurance policy is made inactive, its Default flag is now automatically cleared.
- When studies move from QC to the worklist, they no longer appear in filters that they do not match.
- When selecting to print an image, the print window no longer closes before printing takes place.
- When rescheduling all studies in an order from the reschedule bucket, the appointment confirmation screen now consistently appears.

Summary

- When locally caching a study, priors are now consistently included if Prefetch Priors is set to at least 2 in Rendering Options.
- When Fullscreen is disabled in Viewer Settings, opening a study no longer causes errors.
- When Dragon is selected as the voice recognition option but is not installed, Exa Trans can now be used with Fluency Direct.
- When downloading a 835 payment file, it is no longer downloaded twice, which would result in a disallowed duplicate.
- When DICOM-sending a study, the study is now consistently sent using the transfer syntax configured for the sending AE.
- When certain BTO RCC studies are opened in the viewer, their orientation no longer appears as Medial to Lateral when it should be Head to Foot.
- When adding a study earlier than the number of authorization days, insurance information is no longer erroneously removed from the New Order screen.
- When adding a critical finding from the worklist, a corresponding entry now consistently appears in the audit log.
- When a radiologist adds an additional dictation or addendum to a transcribed or approved study, subsequent users can now consistently see and hear the additional dictation or addendum.
- When a Patient Portal account has multiple authorized representatives, Patient Portal now consistently opens upon sign-in.
- When a combination of single and repeating exclusion schedule rules are set to display as blocks, the user is no longer able to schedule in some of the excluded timeslots.
- Using the Shutter Freehand tool no longer erroneously shuts the entire image.
- Using arrow keys to reposition the Gonstead Ruler no longer causes the ruler to disappear.
- Users are now prevented from directly opening patient details in the patient portal by using a URL without signing in.
- The rendering provider in a claim now consistently defaults to the one configured for the providing facility.
- Routing rules now consistently trigger on study flow from Incomplete to Unread status.
- Querying Modality Worklist no longer fails when the results include a patient whose weight was entered with a decimal value.
- On the viewer toolbar, the appearance of buttons now consistently matches their actual On/Off state.
- Local Cache now skips images with the invalid DICOM tag (0000,0000) rather than attempting to cache them, which caused errors.
- In Web Trans, when editing a new transcription template, the Use Clipboard option is now selected by default.
- In the patient portal, resolved an error occurring when opening the viewer.
- In the Exa Client viewer, when applying C-Spine Analysis annotations, the correct angle now appears without requiring reloading of the viewer.

Summary

- In the Exa Client viewer, restored availability of the Stitch button.
- In the Client Viewer, using the Chiropractic C-Spine Analysis tool no longer causes unresponsiveness.
- In the Claims screen, negative balances are now indicated with parentheses rather than a minus sign.
- In Patient Portal, resolved errors occurring when opening studies in the viewer.
- In Patient Portal, closing an image tab no longer signs the user out of the portal.
- In OpalForwarder, restored the ability to forward files by manually placing them in a drop folder.
- In OpalForwarder, on the Routing tab, destination modalities are no longer case sensitive.
- In Exa Client viewer, when using the CSL tool (from the Chiropractic > Orthospinology menu), the program no longer becomes unresponsive.
- In Billing, when processing an ERA file, all claim adjustments (such as CO 97 and CO 50) are no longer erroneously adjusted off.
- In Billing, the Transaction Detail and Pending Payment screens now consistently display items in chronological order.
- Fixed a problem in which multiple headers appeared when printing images.
- Fixed a problem in which collections processing did not capture all claims that met certain criteria.
- Fixed a problem in which certain events (such as changing critical findings, statuses, flags, or opening exams) were not consistently entered into the audit log.
- An IRS file cache monitoring service no longer consumes excessive system resources when cached studies are updated.
- After searching for a provider on the web, opening the provider record no longer causes errors.
- After merging patient charts, the no longer used chart now appears as inactive on the Recent Patient Results list.
- After installing local cache, it is no longer necessary to sign out and back into the Exa platform to activate it.
- After finding a provider on the web and selecting to add the provider to the Exa platform, updated provider information no longer disappears.
- After DICOM-sending an SR study with an approved report, the name for the created OT file on the worklist is now Approved Report.
- After deleting cached studies and then performing a manual prefetch, the IRS cache indicator (text formatting) on the worklist now consistently updates.